


The hiBBB90 Model Setup Protocol

ver. 1.3

(updated on Mar. 25, 2021)

Overview & Medium composition



Medium name		PGM	PDM ¹	Medium name		AGM	ADM ²
Cells		HBPC/ci37		Cells		HASTR/ci35	
Basal medium		Pericyte medium (Sciencell, #1201-b)		Basal medium		DMEM	
Culture supplements				Culture supplements			
Pen/st (#0503)	1%	+	+	Pen/st	100 units/mL-100 µg/mL	+	+
FBS (#0010)	2%	+	-	N2 (apo)	1%	+	+
Pericyte growth factor (#1252)	1%	+	+	FBS	10%	+	-
Blasticidin S	4 µg/mL	+	-	Blasticidin S	4 µg/mL	+	-
				dBcAMP	1 mM	-	+

¹, Optional. While we used FBS-containing medium in our initial work (Ito et al. 2019), you may change it to FBS-free medium.

², it may be better to further add Glutamax or L-Alanyl-L-Glutamine as we did in our initial work (Ito et al. 2019).


Medium name		BPM
Cells		HASTR/ci35, HBPC/ci37
Basal medium		Neurobasal medium (Life Technologies #21103049)
Culture supplements		
Pen/st	100 units/mL-100 µg/mL	+
N2 (apo)	1%	+
L-glutamine	2 mM	+

Medium name		EM-comp	EM-V/E free
Cells		HBMEC/ci18	
Basal medium		VascuLife BM (#LEB-LM0002) / EBM-2 BM (#CC-3156)	
Culture supplements			
VascuLife VEGF Lfactors (#LEK-LS1020) or EGM-2 SingleQuots Supplements (#CC-4176) ³		+	+
hFGF-b	5 ng/mL	+	+
Ascorbic acid	50 µg/mL	+	+
HC	1 µg/mL	+	+
hIGF-1	15 ng/mL	+	+
hEGF	5 ng/mL	+	-
VEGF	5 ng/mL	+	-
Heparin	0.75 unit/mL	+	+
L-glutamine	10 mM	+	+
FBS	2%	+	+
Pen/st	100 units/mL-100 µg/mL	+	+
Blasticidin S	4 µg/mL	+	-


EM-comp, Endothelial complete medium; EM-V/E free, VEGF- and hEGF-free EM-comp; FBS, fetal bovine serum; HC, hydrocortisone; hFGF-b, human fibroblast growth factor-basis; hIGF-1, human insulin-like growth factor-1; hEGF, human epidermal growth factor; VEGF, vascular endothelial growth factor; Pen/st, penicillin streptomycin; AGM, astrocyte growth medium; ADM, astrocyte differentiation medium; DMEM, Dulbecco's modified Eagle medium; dBcAMP, dibutyryl cyclic adenosine monophosphate; PGM, pericyte growth medium; PDM, pericyte differentiation medium; BPM, brain parenchyma medium.

³, other types of endothelial medium may be used.


Summary of each preparation step (for 12-well format)




Step1-1: Extracellular matrix (ECM) coating




Step1-2: HBPC/ci37 (pericyte) seeding




Step1-3: HASTR/ci35 (astrocyte) seeding




Step2-1: HBPC/ci37 (pericyte) differentiation



Step2-2: HASTR/ci35 (astrocyte) differentiation




Step3-1: HBMEC/ci18 (endothelial cells) seeding




Step3-3: Combine insert and receiver plate

Step3-2: Medium change to BPM



Step1-1: Extracellular matrix coating

Step1-1: Extracellular matrix (ECM) coating



- ❑ Set the inserts on the culture plates^{※1}, then add 250 μ L and 750 μ L ECM solution^{※2} to the apical and the basolateral side, respectively. Stand at 37°C for one hr.
- ❑ Remove the ECM solution and dry the inserts and the plates by air for 30 min.
- ❑ Wash the inserts and the plates with PBS (-).
- ❑ Dry the inserts and the culture plates well for 30 min^{※3}.

※1

In our study, 12-well culture plate and insert (translucent polyethylene terephthalate, 0.4 μ m high-density pores, BD Falcon) are used. However, other culture materials can be used. For 24-well, adjust the volume proportionally.

※2


ECM solution is a mixture of type-IV collagen and fibronectin (final concentration: 0.1 mg/mL of each). Prepare at the time of use.

※3

Make it sure that the inserts are completely dried out.

Step1-2: HBPC/ci37 (pericyte) seeding

Step1-2: HBPC/ci37 (pericyte) seeding



- ❑ Place the insert in a container (which is sufficiently larger than the insert: e.g. 100-mm dish) in an upside-down manner.
- ❑ Prepare 300 μL HBPC/ci37 cell suspension at 3.0×10^4 cells/cm² using the pericyte medium and place them **gently** on the insert^{※4}.
- ❑ Stand it in the clean bench until the cells adhere (approx. 30 min)^{※5}.
- ❑ Turn over the insert upside down and **carefully** set it in the culture plate, where the pericyte medium are filled in advance^{※6}.
- ❑ Culture at **33°C** in 5% CO₂/95% air.

※4

Seed gently, or the cell suspension may fall into the bottom. For a 24-well format, the volume should be adjusted according to the manufacturer's instruction.

※5

No problem if a small portion of cell suspension drops to the bottom during the incubation, unless most cells remain attached.


※6

Confirm that a large number of the cells cannot be found in the bottom of the plate with a microscope.

Step1-3: HASTR/ci35 (astrocyte) seeding


Step1-3: HASTR/ci35 (astrocyte) seeding

Seed HASTR/ci35 cells



Astrocyte growth medium (AGM)

33°C culture for 1 day




- ❑ Prepare 750 μL HASTR/ci35 cell suspension at 2.5×10^4 cells/cm² using the astrocyte growth medium and place it on the culture plate^{※7}.
- ❑ Culture at **33°C** in 5% CO₂/95% air.

※7

Important! To prevent the cells from clustering in the center of the plate, shake gently the plate and make them uniformly attached. For a 24-well format, the volume should be adjusted according to the manufacturer's instruction.

Step2-1: HBPC/ci37 (pericyte) differentiation

Step2-1: HBPC/ci37 (pericyte) differentiation



- ❑ Confirm that a large number of HBPC/ci37 cells are not found in the culture plate with a microscope.
- ❑ Gently wash twice with PBS (-).
- ❑ Change the medium from PGM (for cell proliferation) to PDM (for cell differentiation)^{※8}.
- ❑ Culture at **37°C** in 5% CO₂/95% air^{※9}.

※8


Please note that in our published paper, we use PGM only. While we assume PDM would be better, you can follow the published protocol. In such case, please skip the medium change and just do temperature change alone.

※9

Make sure that the temperature is set at 37°C.

Step2-2: HASTR/ci35 (astrocyte) differentiation

Step2-2: HASTR/ci35 (astrocyte) differentiation



- ❑ Gently wash twice with PBS (-)^{※10}.
- ❑ Change the medium from AGM (for cell proliferation) to ADM (for cell differentiation).
- ❑ Culture at **37°C** in 5% CO₂/95% air^{※11}.
- ❑ Confirm that the HASTR/ci35 cells show differentiated morphology at several hours after the medium change^{※12}.

※10

Important! Wash out FBS completely, or the differentiation would be insufficient.

※11


Make sure that the temperature is set at 37°C.

※12

Please refer to Kitamura et al. J Pharmacol Sci. 2018;137(4):350-358.

Step3-1: HBMEC/ci18 (endothelial cells) seeding

Step3-1:HBMEC/ci18 (endothelial cells) seeding




- Prepare 750 μL HBMEC/ci18 cell suspension at 1.3×10^5 cells/cm² using the VEGF/EGF-free endothelial medium, and gently place it on the insert^{※13}.

※13

For a 24-well format, the volume should be adjusted according to the manufacturer's instruction.

Step3-2: Medium change to BPM

Step3-2:Medium change to BPM




- Change the medium from ADM to BPM (2250 μ L) ※14.

※14

For a 24-well format, the volume should be adjusted according to the manufacturer's instruction.

Step3-3: Tri-coculture

Step3-3:Combine insert and receiver plate



- ❑ Carefully transfer the inserts to the culture plate and start co-culturing.
- ❑ Culture at **33°C** in 5% CO₂/95% air^{※14}.

※14

Make sure that the temperature is set at 33°C.

The hiBBB90 model setup protocol ver. 1.3 Mar. 25, 2021.

This is the first version and will be updated without any prior notices.
Any questions and comments would be greatly appreciated.

Contact Information:

Tomomi Furihata, Ph.D.

Department of Clinical Pharmacy and Experimental Therapeutics

School of Pharmacy

Tokyo University of Pharmacy and Life Sciences

1432-1 Horinouchi, Hachioji, Tokyo, JAPAN 192-0355

Tel/Fax: +81-42-676-8969/5809

tomomif@toyaku.ac.jp